

A2
Cond
of the *IND1* and SEQ ID NO:4 represents the 3' untranslated region of *IND1*. A *IND1* promoter can also be identified by its ability to direct expression in the valve margin of fruit. In particular, the *Ind1* promoter directs expression at the valve margin of developing gynoecium just prior to fertilization (stage 13) through the maturation of the fruit (stage 17). The promoter does not provide significant expression in leaf tissue.--

Please replace the paragraph beginning at page 9, line 29, with the following:

A3
--An "IND1 polypeptide" is a sequence of about 50 to about 200, sometimes 100 to 190, and preferably 198 amino acid residues encoded by a *IND1* polynucleotide. *IND1* polypeptides are characterized by the presence of an basic helix-loop-helix (HLH) domain which bind specific polynucleotide sequences. For instance amino acid residues ISDDPQTVVARRRRERISEKIRILKRIVPGGAKMDTASMLDEAIRYTKFLK (SEQ ID NO:7) represent the HLH domain of the polypeptide shown in SEQ ID NO:2. The HLH domain is known in the art and is shared by other transcription factors including uncharacterized sequences represented by GenBank accession number E1283552 and 2262147 and the gene product, PIF3 (Ni *et al. Cell* 95:657 (1998)). The HLH domain of *IND1* is therefore a DNA binding domain.--

Please replace the paragraph beginning at page 15, line 1, with the following:

A4
--Appropriate primers and probes for identifying genes such as *IND1* from plant tissues are generated from comparisons of the sequences provided herein. For a general overview of PCR see PCR Protocols: A Guide to Methods and Applications. (Innis, M, Gelfand, D., Sninsky, J. and White, T., eds.), Academic Press, San Diego (1990). Appropriate primers for amplification of the genomic region of *IND1* or the *IND1* cDNA include the following primer pairs: 5'-gatgaaaatggaaaatggtatgtata-3' (SEQ ID NO:8) and 5'-gttcatcagggttggagttgtg-3' (SEQ ID NO:9). The amplification conditions are typically as follows. Reaction components: 10 mM Tris-HCl, pH 8.3, 50 mM potassium chloride, 1.5 mM magnesium chloride, 0.001% gelatin, 200 µM dATP, 200 µM dCTP, 200 µM dGTP, 200 µM dTTP, 0.4 µM primers, and 100 units per ml Taq polymerase. Program: 96 C for 3 min., 30 cycles of 96 C for 45 sec., 50 C for 60 sec., 72 for 60 sec, followed by 72 C for 5 min.--